Instructions for Use

OPTIFLUX®

F160NRe, F180NRe, F200NRe, F250NRe
Capillary High Flux Dialyzer
Single Use Only - Electron Beam Sterilized


Figure 1: Dialyzer Orientation During Blood Compartment Priming.

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
Figure 2: Dialyzer Orientation for Hemodialysis Treatment with Arterial-end up.

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
Figure 3: Dialyzer Orientation for Hemodialysis Treatment with Venous-end up.

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GENERAL INFORMATION

Indications for Use: Optiflux Dialyzers are intended for patients with acute or chronic renal failure when conservative therapy is judged to be inadequate.

Contraindications: Optiflux dialyzers are not intended for use in patients that weigh less than 30kg.

USA Only: Federal (US) law restricts this device to sale by or on order of a physician.

Caution: The operator should strictly adhere to the manufacturer's recommended procedures, warnings, and cautions as listed in these instructions for use.

Precautions: Dialyzers may leak resulting in patient blood loss or contamination with dialysate. In the event of a blood leak during dialysis, the health care provider should respond according to the facility's established protocol.

Air entering the extracorporeal circuit during dialysis can result in serious injury or death. Check the security of all extracorporeal connections prior to the initiation of dialysis and periodically throughout the treatment. The venous drip chamber should be continuously monitored with a level / air detector.

Warning: Due to the high water flux capability of high permeability membranes with an ultrafiltration coefficient > 8, it is necessary to use such dialyzers only in conjunction with dialysis machines that are equipped with precise ultrafiltration control, such as the Fresenius 2008 series. In any case, the safety instructions for the hemodialysis machine must be followed.

Dialyzer performance may change with the duration of observation.

The user is cautioned to regularly monitor the patient's chemistry values using quantitative measurements and analysis to ensure that the prescribed therapy is delivered. The clinical parameters monitored should at least include urea, hemoglobin and serum albumin.

Dialysate: The dialysate must meet current water quality standards for dialysis (ISO 11663) to mitigate the risk of infection due to back-filtration.

Side effects: Adverse reactions such as hypertension, hypotension, headache, and nausea may be associated with hypovolemia or hypervolemia and can usually be avoided with careful management of patient fluid, electrolyte balance, blood flow rate, and ultrafiltration rate.

In rare cases, thrombocytopenia or hypersensitivity reactions including anaphylactic or anaphylactoid reactions to the dialyzer, or other elements in the extracorporeal circuit, may occur during hemodialysis. Hypersensitivity reactions may cause mild to severe signs and symptoms, including: itching, flushing, hives, swelling, fever, leukopenia, hypotension, hypertension, shortness of breath with wheezing, arrhythmias, and/or respiratory arrest.

Patients with a history of hypersensitivity reactions or patients who have a history of being highly sensitive and allergic to a variety of substances should be carefully monitored during treatment. With severe hypersensitivity reactions, dialysis must be discontinued and aggressive first line therapy for hypersensitivity reactions must be initiated. The decision to return the patient's blood in the event of a hypersensitivity reaction is determined by the physician.

Heparinization: We recommend systemic heparinization of the patient. Systemic heparinization is the administration of the prescribed loading dose of heparin into the vascular access of the patient; then waiting a period of 3 to 5 minutes prior to initiating the treatment. During dialysis, the dose of heparin and method of administration is the decision of the physician.

Sterile/Non-pyrogenic: The blood pathway is sterile and non-pyrogenic if blood port caps are in place and undamaged. Do not use if the dialyzer is damaged in any way. Use aseptic technique for all blood side connections. Structural integrity of the hemodialyzer is warranted for the first use only when prepared as directed. **Sterilized using the electron beam method of sterilization.**

Recommended storage: Between 5 and 27 degrees C (41 - 80 degrees F).

Dialyzer reuse: Optiflux dialyzers are not designed for or intended for reuse.

Note: Follow the hemodialysis machine specific operating instructions for treatment, when applicable; otherwise, refer to the instructions below.

PREPARATION FOR DIALYSIS

• If the dialysate delivery system was chemically disinfected or sterilized prior to patient use, be sure to test for the absence of germicide residuals with a test intended for this application, according to the test manufacturer's instructions.

• Place the dialyzer in the dialyzer holder in the vertical position.

• Install the arterial and venous bloodlines on the hemodialysis machine.

Note: Refer to dialysate delivery machine manufacturer's instructions for use for setting up bloodlines.

• **Remove the dialyzer cap from the arterial end of the dialyzer and immediately secure it to the dialysate port. Aseptically connect the arterial (red) dialyzer end of the bloodline to the lower dialyzer blood port.**

• **Remove the dialyzer cap from the venous end of the dialyzer and immediately secure it to the dialysate port. Aseptically connect the venous (blue) dialyzer end of the bloodline to the upper dialyzer blood port.**

• Ensure all connections are secure. Refer to Figure 1: Position of Dialyzer During Blood Compartment Priming.

• Aseptically spike a 1 liter bag of 0.9% sterile saline solution with a clamped dialysis priming set.

• If not already attached, attach the dialysis priming set to the saline "T" connection located just before the blood pump segment on the arterial bloodline. Check to be sure the connection is secure.

• Open the clamp on the dialysis priming set and allow saline to gravity prime the portion of the arterial bloodline from the saline "T" to the patient end.

• Clamp the main line tubing on the arterial bloodline between the patient end and the saline "T" connection.

• Start the blood pump and set a pump speed of 150mL/min. Prime the remainder of the arterial bloodline, dialyzer, and venous bloodline with saline. While the extracorporeal circuit is filling with saline, intermittently pinch and release the bloodline between the blood pump and the dialyzer to help purge air from the dialyzer. Gently tap the dialyzer to facilitate air removal from the dialyzer.

• Fill the dialyzer and bloodlines with 300mL sterile 0.9% saline solution. The drip chambers in the bloodlines should be set to and maintained at 3/4 full.

• Stop the blood pump. Clamp the venous bloodline. Aseptically connect the patient ends of the arterial and venous bloodlines together in preparation for recirculation of the extracorporeal circuit. Unclamp main line clamps on arterial and venous bloodlines.

• Perform Pressure and Alarm tests.

• Verify that the dialysate is within the prescribed conductivity limits with a calibrated independent conductivity meter. To identify situations where the acetate or acid and bicarbonate concentrates are not properly matched, use a calibrated pH meter to verify that the pH of the dialysate is within the appropriate physiologic range.

• Remove and save the dialyzer caps, and attach the dialysate lines to the dialysate ports of the dialyzer. In order to maximize the efficiency of the dialyzer assure that the dialysate lines are attached so that the blood and dialysate flows are counter-current.

• Rotate the dialyzer so the venous end is lower than the arterial end to fill the dialysate compartment with dialysate.

• Assure that the dialysate compartment is completely filled with dialysate.

• Return the dialyzer to the arterial end down position.

• Recirculate the extracorporeal circuit at a blood flow rate of 300 to 400mL/min and a dialysate flow 500mL/min until all air has been purged from the dialyzer and bloodlines.

• During recirculation, to assist in removing air from the dialyzer, intermittently pinch and release the blood tubing between the blood pump and the dialyzer. It is recommended to gently tap the upper portion of the dialyzer to remove air even if the header area of the dialyzer looks free of air.

• Do not infuse the recirculated saline prime into the patient. Discard the recirculated saline and fill the entire extracorporeal circuit with fresh saline prior to connecting to the patient. The volume of fresh saline used to fill the extracorporeal circuit should be equal to the volume of the dialyzer and blood tubing set in use.

INITIATION OF DIALYSIS

• To initiate dialysis; reduce the blood pump speed to 150 mL/min and then stop the blood pump, clamp the dialysis priming set, arterial, and venous bloodlines.

• Aseptically attach the patient ends of the bloodlines to the arterial and venous access of the patient. Open the arterial and venous mainline and vascular access clamps.

• Turn on the blood pump and slowly increase the blood pump speed to the prescribed blood flow rate. Be sure to monitor the hemodialysis machine arterial and venous blood pressures carefully during this process to note any possible flow restrictions or inappropriate pressure readings.

• Once the prescribed blood flow rate has been achieved, set the prescribed ultrafiltration rate.

• Performance testing has demonstrated that the dialyzer can be run in either the arterial end up or arterial end down position during the hemodialysis treatment. Refer to Figure 2 or 3: Hemodialysis Treatment for proper dialyzer and bloodline positioning.

DURING THE DIALYSIS TREATMENT

• If a blood leak should occur during the treatment, the operator should follow established facility procedures for a dialyzer blood leak.

• Air entering the extracorporeal circuit during dialysis is a very serious event and should be avoided. A routine check of all connections prior to initiation of dialysis and periodically throughout the dialysis treatment is recommended. Constant monitoring of the venous drip chamber with a level / air detector is required. Should air enter the venous line during the treatment, the dialysis treatment must be discontinued without returning any of the blood mixed with air.

Instructions for Use

OPTIFLUX®

F160NRe, F180NRe, F200NRe, F250NRe
Capillary High Flux Dialyzer
Single Use Only - Electron Beam Sterilized

COMPLETION OF DIALYSIS

• When the dialysis treatment is completed, turn the blood pump off and set the UF rate to the recommended minimum. Verify sufficient 0.9% sterile saline solution for return of the blood in the extracorporeal circuit back to the patient.

• Open the priming line clamp and place a hemostat on the arterial bloodline between the saline "T" and the blood pump. Return the blood in the mainline between the saline "T" and the patient end back to the patient.

• Close the arterial mainline clamp when return is completed.

• Remove the hemostat and start the blood pump at 100-150 mL/min pump speed. Intermittently pinch and release the blood tubing between the blood pump and the dialyzer to help to efficiently return the blood in the extracorporeal circuit back to the patient. Do not let air enter the extracorporeal circuit.

• Once the blood has been returned to the patient, turn the blood pump off. Clamp the venous mainline and the patient's arterial and venous access. Aseptically disconnect the arterial and venous bloodlines from the patient's access.

• Remove the dialysate lines from the dialyzer and apply the dialyzer caps.

Note: If dialyzer caps are not available, rotate the dialyzer horizontally making sure to have the dialysate ports facing up to prevent leaking before removing the dialysate lines from the dialyzer.

• Discard the extracorporeal circuit in an appropriate biohazard waste receptacle. References: 29CFR, 1910.145, 1910.1030 (Code of Federal Regulations) and appropriate state and local codes.

Technical Data: These data represent typical in vitro performance. Actual in vivo performance may differ.

	Optiflux F160NRe	Optiflux F180NRe	Optiflux F200NRe	Optiflux F250NRe
Catalogue Number	0500316E	0500318E	0500320E	0500325E
Membrane Surface Area (m²)	1.5	1.7	1.9	2.5
Ultrafiltration Coefficient¹ (mL/hr/mmHg), Q _b =300mL/min	61	76	74	111
Priming Volume Blood (mL)	87	102	113	142
K_oA (mL/min) Q _b =300mL/min, Q _d =500 mL/min	1167	1321	1415	1714
Flow Resistance (mmHg)				
Blood Q_b=300mL/min	89	80	66	57
Dialysate Q_d=500mL/min	17	15	15	19
Maximum TMP (mmHg)	600	600	600	600
Minimum Blood Flow (mL/min)	100	100	100	100
Maximum Blood Flow (mL/min)	600	600	600	600
Minimum Dialysate Flow (mL/min)	200	200	200	200
Maximum Dialysate Flow (mL/min)	800	800	800	800
Maximum Ultrafiltration Rate (mL/hr)	4000	4000	4000	4000
Clearance² (mL/min) Q _b =300mL/min, Q _d =500mL/min,				
Urea ³	271	277	280	287
Creatinine	242	253	259	273
Phosphate	241	253	254	274
Vitamin B ₁₂	153	169	177	205
Lysozyme ⁴	74	87	83	113
Q_b=400mL/min, Q_d=500mL/min,				
Urea ³	317	328	333	349
Creatinine	274	290	298	321
Phosphate	279	297	298	328
Vitamin B ₁₂	162	181	189	226
Q_b=400mL/min, Q_d=600mL/min,				
Urea ³	332	344	348	366
Creatinine	288	304	310	338
Phosphate	290	307	309	341
Vitamin B ₁₂	167	189	195	235
Q_b=400mL/min, Q_d=800mL/min,				
Urea ³	349	360	365	381
Creatinine	302	320	328	357
Phosphate	304	321	324	356
Vitamin B ₁₂	173	195	203	248
Sieving Coefficient⁵ Q _b =300mL/min, ultrafiltration rate=29mL/min				
Albumin			<0.01	

Dialyzer Technical Information

Membrane Material:	Advanced Fresenius Polysulfone®
Fiber Inner Diameter:	185 microns
Membrane Wall Thickness -- F160NRe, F180NRe, F200NRe:	35 microns
Membrane Wall Thickness -- F250NRe:	30 microns
Housing Material:	Polycarbonate
Potting Compound:	Polyurethane
O-ring:	Silicone
Recommended Blood Connections:	ISO 8637-1:2017
Recommended Dialysis Fluid Connections:	ISO 8637-1:2017
Sterilization Method:	Electron Beam

Footnotes:

1 Bovine blood, Hct 32%, protein 60g/L, 37°C

2 Aqueous solution of respective test substance, 37°C; ultrafiltration rate=0mL/min

3 Sodium used as a marker for urea

4 Lysozyme used as a marker for middle molecule solutes; molecular weight 14,300 Daltons

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